Instruction for use
HEINE DELTA[®]one
HEINE DELTA[®] 30
HEINE DELTA[®] 30 PRO

HEINE DELTA[®]one

HEINE DELTA[®] 30

HEINE DELTA[®] 30 PRO

DEUTSCH

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Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.	
	The supply voltage quality should be that of a typical hospital environment.	
	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.	
Performance features of the ME system that have been determined to be essential to the performance	None	
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle		
Warning	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.	
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.	
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the considered device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.	
	DELTA 30/DELTA 30 PRO/DELTAone: Illumination failures may occur due to electrostatic discharge.	
Note	The EMISSIONS characteristics of the DELTA 30 PRO make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), the DELTA 30 PRO might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the DELTA 30 PRO.	
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used as described in the chapter "Options". The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.	
Test	Compliance	
RF emissions CISPR11	DELTA 30 / DELTAone	DELTA 30 PRO
	Group 1 Class B	Group 1 Class A
Conducted emissions (EN 55011/CISPR 11)*	Passed	
Radiated emissions (EN 55011/CISPR 11)		
Harmonic current emissions (IEC 61000-3-2)*		
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*		
Immunity	See attached immunity test levels	

Immunity test levels

Test	Test level	
	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ± 8 kV Air Discharge: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m 80 – 2700 MHz 80 % AM at 1 kHz	
Electrical fast transients / bursts (IEC 61000-4-4)*	± 2 kV 100 kHz repetition frequency	
Surges (IEC 61000-4-5)*	± 0.5 kV, ± 1 kV	
Conducted disturbances induced by RF fields (IEC 61000-4-6)*	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM: ± 5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 50 Hz or 60 Hz	DELTA 30/DELTA 30 PRO: Compliance test levels are identically to the IEC 60601-1-2 test levels DELTAone: Test was not necessary, because there are no magnetic sensitive components.
Voltage dips (IEC 61000-4-11)*	0 % U_T ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles Single phase: at 0°	DELTA 30/DELTA 30 PRO: Compliance test levels are identically to the IEC 60601-1-2 test levels DELTAone: Test was not performed, because it was not considered necessary. No risks were identified.
Short interruptions (IEC 61000-4-11)*	0 % U_T ; 250/300 cycles	

*n/a: "Not applicable" in the internally powered mode

Technical specification

Environmental specification DELTAone, DELTA 30 and DELTA 30 PRO

Environmental conditions for operation	DELTA 30/DELTAone +10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa	DELTA 30 PRO +10 °C to +35 °C 10 % to 75 % non-condensing 700 hPa to 1060 hPa
Environmental conditions for storage	DELTA 30/DELTAone +5 °C to +45 °C 45 % to 80 % rel. humidity	DELTA 30 PRO Less than 1 month: +5 °C to +50 °C Less than 3 months: +5 °C to +40 °C Less than 1 year: +5 °C to +20 °C rel. humidity < 70 % 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa	

Technical specification DELTAone

Power supply	USB 2.0 Micro B: 5 V DC, 650 mA
Power consumption	max. 3.25 W
Class	Charging: II Operation: internally powered
Charging time	max. 75 min
Operating time	typ. 80 min at highest brightness in polarisation mode
Device classification according to IEC 62471	Risk group 1
Applied part	Contact plate type BF

Technical specification DELTA 30 and DELTA 30 PRO


Power supply	DELTA 30 USB 2.0 Micro B: 5 V DC, 1.2 A HEINE Charger 30	DELTA 30 PRO USB-C: 5 V DC, 1.2 A HEINE Charger 30
Power consumption	max. 6 W	max. 6 W
Class	Charging: II Operation: internally powered	Charging: II Operation: internally powered
Charging time	USB: typ. 120 min. HEINE Charger 30 typ. 120 min.	USB-C: typ. 150 min. HEINE Charger 30 typ. 150 min
Operating time	typ. 160 min at highest brightness in polarisation mode	Min. 210 min. @ 5300 K, Polarised, 100 % Brightness
Device classification according to IEC 62471	Risk group 1	Risk group 1
Applied part	Contact plate, working ring type BF	Contact plate, working ring type BF

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HEINE Optotechnik GmbH & Co. KG