
Hikoneb 908 DC

User Manuel EN

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- Large Ultrasonic Nebulizer
- Ultraschallvernebler

CE 1984

Rev. 03_2023

LARGE ULTRASONIC NEBULIZER

Table of contents

1. Symbol Definitions	14
2. Important Safeguards	15-17
3. How the device works	17-18
4. Device Set-Up	18-20
5. Device Image and Components	20-21
6. Operating the device	22-23
7. Device Settings	23-25
8. Troubleshooting	26
9. Maintenance, Repair and Disposal	27-28
10. Cleaning, Disinfection, Sterilization	28-30
11. Accessories	30
12. Technical Specifications	31
13. Electromagnetic Conformity	32-34

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1. Symbol Definitions

	It is mandatory to read and understand the operating instructions prior to use. ① This symbol has a blue background on the product label		This device contains electronical and/or electrical equipment that must be recycled per EU Directive 2012/19 EU-Waste Electrical and Electronic Equipment (WEEE)
	Type of Protection. Ground Conductor as per Class I		Medical device
	Keep dry		IAC Voltage
IP20	The device is protected against solid objects of 12.5mm and more		Item code
	Manufacturer		Date of production
	Type B applied part		European CE-Mark
	Serial Number		

2. Important Safeguards

When using electrical products, especially when children are present, basic safety precautions should always be followed.

Read all instructions before using. Important information is highlighted by these terms:

DANGER – Urgent safety information for hazards that will cause serious injury or death.

WARNING – Important safety information for hazards that might cause serious injury.

CAUTION – Information for preventing damage to the product.

NOTE – Information to which you should pay special attention.

DANGER!

To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a bath tub or sink.
4. Do not place in or drop into water or other liquid.

WARNING!

To reduce the risk of burns, electrocution, fire, or injury to persons:

1. The device should never be left unattended when plugged in. Close supervision is necessary when this product is used for or on near children or physically incapacitated individuals.
2. Use this product only for its intended use as described in this guide. Only use attachments/spare parts recommended by the manufacturer.
3. Do not nebulize home-therapy or aroma-therapy oils with this device.
4. Do not use on ventilated patients.

Never operate this product if:

1. It has a damaged power cord or plug.
2. It is not working properly.
3. It has been dropped or damaged.
4. It has been dropped into water.

Return the product to an authorized Drive DeVilbiss service centre for a technical inspection and repair.

WARNING!

- Keep the power cord away from heated surfaces.
- Never drop or insert any object into any opening or tube.
- Do not use outdoors. This device is designed to be used inside.
- Do not use in an oxygen enriched environment.
- Connect this product to an appropriate electrical outlet only.
- Do not use this product with a DC to AC converter, or with any AC voltage and/or frequency other than specified.
- The Hikoneb 908 complies with EMC requirements conforming to IEC:60601-1-2. Nevertheless, location should be considered with caution, when this product is used in the vicinity of other electronic equipment where electromagnetic interference may be an issue.
- Make sure that all parts are connected correctly, and that the dust filter is clean.
- A filter which has changed its colour or used a period of 3 months, should be replaced.
- Bacteria filter should be replaced every 48 hours, in between patients and when it gets dirty.

NOTE!

- Medicine cups (38) are for single use only.
- Dispose of any remaining medicines. Do not use again.
- The device has been designed for inhalation using a mouthpiece or a face mask. When only using the aerosol tube (for humidification purposes), there should be at least 10cm distance between a patient and the tube.
- Change mouthpiece daily or when changing patient.
- Room/ambient temperature should not exceed 30°C when used with the heater system (19).
- Only medical clean/sterile water should be used for nebulization.
- Medicines can be nebulized as aqueous solution only.
- The nebulizer chamber (11) is equipped with a maximum and minimum water level mark. Water should be filled into the chamber between these two marks. If a sterile water bottle is used, the gasket should be placed completely on the bottle lid, to prevent flowing water from the chamber. In addition, make sure that the lid is tightly

secured on the sterile water bottle.

- After each use, this system should be cleaned, disinfected, and sterilized.
- The device should be used with original manufacturers air filter (23).
- Caution should be taken when removing the chamber, tubing/tube heater after therapy, as these items may be very warm.
- Only authorized service personnel are allowed to open the device to conduct troubleshooting and failure analysis. Do not try to open or repair the device by yourself.
- Use original parts and accessories only.

3. How the device works

The Hikoneb 908 uses an ultrasonic frequency to turn a liquid medication into an aerosol.

The device can be started and controlled by using push buttons. The base of the nebulizer chamber contains a transducer, which causes the liquid in the chamber to nebulize through its high-frequency oscillations. A blower located in the housing conveys the aerosol through the patient tube to the patient for inhalation. This driving air is passed through a bacteria filter before nebulization.

The Hikoneb 908 is equipped with an alarm system. An audible alarm is activated when water in the nebulization chamber is below the minimum water level and the device will automatically switch off.

Service Life

The expected service life of the device is 10 years.

Intended Use / Indication

The Hikoneb 908 DC has been manufactured for Ultrasonic Nebulizer Aerosol treatment. It breaks water or liquid medicine into particles and produced respirable water or medicine steam. It is used for Asthma, Bronchitis, Pharyngitis and Allergies and other respiratory diseases.

The Hikoneb 908 DC is used in hospital environments, hospital employees make usage and maintenance operation. It is not suitable to connect to anesthesia respiration system or ventilator systems. A doctor decides medicine type and

usage, a user cannot practice this treatment directly by himself/herself. The devices are not support, patient monitoring and sterilized devices. The device has no measurement function.

Even though high-volume ultrasonic nebulizers are not intended for direct patient application, mouthpieces and aerosol masks can be used when administering medications and substances intended to be inhaled (antibiotics, broncho-dilators and for sputum induction)

Not suitable to be used with respiration masks. This warning is to avoid using with protective masks such as N95. If a mask is to be used, then use aerosol/nebulizer masks only.

Device intended use is for humidification of patient's respiratory air. Room humidification for medical and non-medical applications use is prohibited. Manufacturer is not liable for product failure, all kind of damages and accidents caused by non-intended use.

Contraindications

- Coughing, bronchospasm or allergic reactions.
- Patients unable to spontaneously breathe.
- Inability of the patient to cooperate accordingly.

Target Group

Spontaneously breathing patients from the age of 16 years for aerosol therapy.

4. Device Set-up

Assembly Instructions

- 1.** Assemble the 5 feet mobile stand by using a 6mm Allen key (included) to connect the five wheels and five feet. (26 & 1) with the black middle base connector (3) and 6mm Allen bolts. Screw together the main stand rods (6) and connect to the base with a 6mm bolt and Allen key.
- 2.** Connect the main unit bracket clamp (5) with the assembled stand rod (6) by its clamp and Allen key. Connect three knurled plastic locking nuts (4) below the device,


to the device's (5). The distance between the device and the ground should not exceed 60cm; a higher position can affect stability of the assembly.

3. Clip and fix the bottle holder (8) to the stand (6) with its screw.

4. If entraining oxygen, connect the oxygen adapter to device outlet, then connect the bacteria filter (23) to the oxygen adapter (21) and connect the incoming oxygen to O2 adapter. Oxygen adapter must always be connected before the bacterial filter (between device outlet and bacterial filter). If connected between bacterial filter and ventilation tube, oxygen will enter to the system without any bacterial filtration.

5. If heated nebulization is desired, then place the heater (19) on chamber between the lid and the aerosol tubing (15) and connect its power plug to the socket on the back of the device.

BACK SIDE OF THE DEVICE


6. Attach one end of the ventilation tube (16) to the bacteria filter (23) and the other end to the nebulization chamber lid.
7. Connect the adjustable arm set (14) to the main stand rod. Connect the aerosol tube (15) to the arm set's black plastic clips. The support arm should tighten and secured ensuring that the position of the aerosol tube is pointing upwards and straight, allowing any collected aerosol to flow back into the chamber.
8. To start the aerosol therapy, place the required dose of medication into the disposable medication cup (38), and connect to the chamber (11) and its lid (36). Locate and lock the nebulization chamber in its place to the unit (22) by turning to the right.

5. Device Image and Components

1. Wheel
3. Middle base connector
4. Locking nuts
5. Bracket clamp
6. Stand rod
8. Bottle holder
11. Nebuliser chamber
14. Adjustable arm set
15. Aerosol tube
16. Ventilation tube
19. Heater
21. Oxygen adapter
22. Holding device
23. Bacterial filter (round)
26. Five feet
36. Nebuliser chamber lid
38. Medical cup

NOTE—The distance between device and ground should not exceed 60 cm.

6. Operating the device

Nebulization directly from the nebuliser chamber

After the device is assembled use a sterilized nebulizer system and ensure that both, air filter (33) and bacteria filter (23), are properly installed.


- Medication in doses between 5ml and 10ml can be nebulized.
- Remove the lid (36) from the nebulization chamber (11) by turning it slightly back and forth and pulling at the same time.
- Plug the conical medication cup (38) on the nebulization chamber lid (36) make sure the water is between the minimum and maximum levels.
- Place the nebulizing chamber lid to the nebulizing chamber (11). Insert the nebuliser chamber into the holding device (22) and lock the nebuliser chamber by turning it to the left until it clicks into place.
- Apply desired medication (max 10ml) into the disposable medication cup (38).
- Connect the aerosol tube (15) and ventilation tube (16) to the chamber lid.
- Turn on the device using the ON/OFF switch which is located on the left side of the control panel.
- Switch on the device with the mains switch (I). The nebulization process will start.

NOTE!

If the nebulization process does not start immediately, check that the nebulization chamber is correctly seated in the device holder.

NOTE!


Do not let the nebulizer chamber run dry. Low Level Water alarm: If the water in the Nebulization Chamber falls below the lower limit, the warning led will illuminate and the unit will alarm. When water is added, both warning LED and alarm will deactivate.

7. Device Settings

The Hikoneb 908 ultrasonic nebuliser has various menu settings in order to accommodate individual patient requirements. These menu settings consist of a Timer Setting, Aerosol Volume, Ventilation Speed and a Heater function. All settings can be accessed via the "Menu"-button. A green LED indicates the position within the Menu..

Operation Panel

Timer Setting

By default, the Timer Setting is set to continuous operation C. By pushing the + button the integrated timer can be set between 1 and 9, which is equal to 10 – 90 minutes. 0 stops the operation of the device, all other settings will remain. Once the Timer Setting has been activated, the timer will count down (in steps of 10) to 0.

Aerosol Volume

By default, the Aerosol volume is set to 2. The Aerosol Volume density is adjustable in steps between 1 and 4. By pushing the – or + button the Aerosol Volume produced, can be increased, or decreased.

Ventilation Speed

By default, the Ventilation Speed is set to 2. The Ventilation Speed is adjustable in steps between 1 and 4. By pushing the – or + button the nebulization flow can be increased or decreased.

Heater Function

By default, the Heater Function is set to 0. Pushing the + button will activate the heater. An on/off LED (I/O) is illuminated when set, to let the user know that the Heater Function has been activated.

NOTE!

- Before practicing medicines, product information provided by the manufacturer should be examined carefully. A qualified physician should administer/prescribe the dosage of medicine to be used.
- Only liquid medicine/solutions shall be used.
- Maximum medication 10 ml.
- Therapy Process and respiration period should be determined by a qualified physician.
- Due to its chemical characteristics, the nebulization chamber, medication cup and aerosol & ventilation tubes are not suitable for ether type lubricants, ketones, hydrocarbons, and cyclohexanone type substances.

8. Troubleshooting

If a fault occurs, first switch off the device. To do this, press the power button. If possible, eliminate the cause of the fault and switch the device on again.

Error	Reason	Solution
The device is not operational, On-Off button light is not on	No mains	Check plug, IEC connections and fuses
On-Off button is on but there is no aerosol production	1. Transducer 2. O/C fuse	1. Replace transducer 2. Replace fuses
Red error LED is on	Low Water Level Nebulizing chamber not locked on the device correctly	1. Add more water 2. Re-fit chamber ensuring correct
Aerosol production is visible, but not coming out of tubing	Blockage of patient tube due to wrong positioning Occluded Filter Flow menu disabled Faulty fan	Change position of the aerosol tube to point upwards and straight Replace the filter Activate Flow setting Replace Fan

8a. Replacement of the main fuse.

WARNING!

Electric shock! Before replacing the fuses, disconnect the device from the mains.

The main fuse (2A) is located inside the IEC socket on the back of the device. Remove the power cord from the device. To remove fuse holder, use a small flat screwdriver pulling the holder forward. Replace fuse with the same fuse rating and replace fuse holder by pushing back into position.

9. Maintenance, Repair and Disposal

Looking after your device and keeping it in good condition will ensure the device lasts the expected service life. To keep the device in good working order, we advise to clean the device after every use, and would suggest functional testing every six months as per the instructions detailed in this user manual. Should the device ever become faulty, technical support may be required.

Bacteria filters can be disposed of with normal domestic waste. Nebulization chamber, mobile trolley wheels, aerosol tubing which include PC, ABS and PVC can also be considered as normal waste. The outer casing of the device itself is constructed from aluminium and would need to be recycled or disposed of under WEEE regulations.


Devices returned to the manufacturer for disposal, will be processed under WEEE Regulations.

Bacteria Filter Replacement

Using the approved bacteria filter (23) is mandatory to prevent contamination during aerosol therapy. The filter should be replaced after 48 hours of continuous use, or 1 week for non-continuous use. Low aerosol output is also a good indication of a contaminated/occluded filter.

Transducer Replacement

Turn off the device, remove the Nebulization chamber (9) and take off the lid of nebulization chamber. Undo the chromium-plated nut underneath the nebulization chamber and remove the complete transducer assembly (41). Replace the new transducer so that the shiny side is facing towards the nebulization chamber. Reapply and tighten the chromium nut.

**NOTE!**

The device will not operate if the transducer is not installed correctly. Also, whenever the transducer has been replaced and before the chamber is connected in place, fill the chamber with water and check for possible leaks.

10. Cleaning, Disinfection, Sterilization

Cleaning

The device and 5 feet mobile stand can be cleaned with a mild detergent and damp cloth. For special cleaning/disinfecting please follow the manufacturers cleaning instructions detailed in this user manual.

NO parts should be exposed to excessive pressure during the cleaning process. Pay specific attention to the transducer surface ensuring it is clean. Clean lime scale residue by wiping, not scrubbing. Abrasive scrubbing can damage the transducer.

Disinfection

When disinfecting the device surface, the device's surface, make sure that the mains power cable is removed. Ensure that no liquid enters inside the device when cleaning. Operate the device only after the disinfection has completely dried and evaporated. Use only disinfection products consisting of aldehyde, ammonia and alcohol components. Never use phenol derivative disinfections, as this will compromise and reduce the service life of plastic parts.

Examples of surface disinfectants:

- Incidin® Plus *
- Incidin® Perfect *
- Antiseptica Kombi surface disinfectant *
- Antifect® FF *
- B 10 *

Examples of surface disinfectants:

- Cidex® AF
- Gigasept® FF *

* Registered trademarks of the respective manufacturers

The lifespan of the nebulizer chamber with lid and the quartz module is limited to 30 autoclaving cycles. The example disinfectants are recommendations from the device manufacturer. The recommendation of the mentioned active ingredient bases (pH 4.5-8) applies. This does not release the user from the need to consult a hygiene specialist in individual cases or to test the components to be disinfected for compatibility with the active ingredients.

Sterilization

Certain parts of the device and nebulizer chamber can be steam sterilized. Steam sterilization 134°C one cycle, 3 minutes duration is allowed only for the following parts:

Nebulizer Chamber (11)
Nebulizer Chamber lid (36)
Aerosol Tube (15)
Ventilation Tube (16)

Before sterilizing, ensure that any cleaning or disinfection residue is completely removed. Once sterilized, please examine all parts for possible damage. Any damaged part should be replaced before use.

NOTE!

Disposable Medication Cup (38) is for single use only.

WARNING!

Risk of injury! Disinfectants can contain harmful substances that can cause injuries when they come into contact with skin and eyes. Protect your skin and eyes and observe the rules of hygiene when working with disinfectants.

11. Accessories

Item code	Description	Figure
X000658	Air filter (pack of 2)	33
M000042	Nebulizer Chamber, complete (includes lid and transducer)	11, 36, 39
M000436	1,7Mhz transducer	39
X000661	Bacterial filter (round)	23
X000659	Ventilation tube	16
X000660	Aerosol tube	15
X000634	Oxygen adapter	21
M000046	Heater	19
G-252-001	Power plug (3 prongs), EU	
M000038	Medical cup, disposable (x50)	38
X002496	Bottle hanger (plastic)	8
M000061	Set of 5 feet mobile stand and support arm	1, 3, 6

12. Technical Specifications

Class	I
Type	B
Power	230V, 50-60Hz, 0.21A
Transducer Frequency	1.7 MHz
Dimensions	290 x 226 x 108mm
Weight	3kg
IP	20
Operation Conditions	5/40°C temperature, 10%/95 non-condensing
Transport Conditions	-30/70°C temperature, 15%/95 non-condensing
Operation Pressure and Flow	>60Kpa, >4L/dk
Sound Pressure Level	<65dB
MMAD % particles < 5 µm	4.22 µm 60%
Nebulization Speed	2g/min.
Aerosol Output Rate	max. 2.3ml/min.
Operating mode	continuous operation, timer
Classification according to EC directive 93/42/EEC Annex IX	Ila

Nebulization Chamber

Weight	0.407kg
Internal Volume	0.743L
Operation Volume	0.277 – 0.348L

13. Electromagnetic Conformity

The nebulizer can operate without being affected by electromagnetic frequencies at a certain level. Ultrasonic Nebulizers Electromagnetic Conformity Directive (89/336/AT). The device has been tested and manufactured conforming to the Medical Devices EEC/93/42 Medical Devices Rules and EN 60601-1 and EN60601-1-2 Electromagnetic conformity. EMC Test

Test	Level	Ambient
TS EN 61000-4-2 Electrostatic Discharge Immunity	With contact: ± 8 kV From air: $\pm 2 - 4 - 8 - 15$ kV	Laying where the device is placed, should be wooden, concrete or ceramic. If laying is synthetic, relative humidity should be 30% at least.
TS EN 61000-4-3 Beaming, RF, Electro Magnetic Area Immunity	10V	Portable and movable RF devices should not be approached to our device. Required distances have been given below.
TS EN 61000-4-4 Electrical Rapid Temporary Regime/ Explosion Immunity	± 2 kV 100 kHz	Power feeding should be on mains or hospital installation quality, in building power input, Leakage Current Protective Relay should be plugged if there is not.
TS EN 61000-4-5 Immunity Against Sudden Increases	L-N: $\pm 0,5 - 1 - 2$ kV L-PE: $\pm 0,5 - 1 - 2$ kV	Power feeding should be on mains or hospital installation quality, in building power input, Leakage Current Protective Relay should be plugged if there is not.
TS EN 61000-4-6 Its Immunity Against Inducted and Transmitted Breaks by RF Areas	3 V	Portable and movable RF devices should not be approached to our device. Required distances have been given below.
TS EN 61000-4-8 Main frequency magnetic area immunity	30 mA	Its power frequency magnetic area should be in hospital or mains quality. Our devices have no magnetic break effect on power network.

Test	Level	Ambient
TS EN 61000-4-11 Voltage Dips, Short Interruptions and Voltage Changes Immunity	%0- 0,5 time 0°-45°-90°-135°-180°- 225°270°-315° %0 – 1 time %70- 25/30 times 0° single phase	Power feeding should be on mains or hospital installation quality, in building power input, Leakage Current Protective Relay should be plugged if there is not. If a user wants to use the device in electric cuts, it is required to use Uninterruptable Power Supply.
TS EN 55011 Terminals disruptive Voltage Emission Break	Class B	RF transmission is at low level, it does not jam in other devices.
TS EN 61000-3-2 Harmonics	Class A	It is suitable to be operated at Low Voltages in mains and hospital installation type, it can directly be connected to a socket. If Instantaneous Voltage changes occur often, Voltage Regulator is recommended.
TS EN 61000-3-3 Voltage Fluctuations and Flicking		

Distances recommended for preventing RF frequencies

Placing the devices within the following distances from fixed and portable RF devices will assist in preventing interferences.

Nominal power of the transmitter	150KHz-80MHz ISM bands and 80MHz-800MHz $u=0.5\sqrt{G}$	800MHz-2.5GHz $u=2.3\sqrt{G}$
0.01 W	0.12 m	0.24 m
0.1 W	0.37 m	0.74 m
1 W	1.17 m	2.34 m
10 W	3.69 m	7.38 m
100 W	11.67 m	23.34 m

u = distance– meter, G = Power– Watt

ISM band: Industrial, Scientific, and Medical bands.

Please refer to the table above when choosing a location and safe distance from RF devices to operate the ultrasonic.

These distances might not be valid in all situations. Electromagnetic spread is also affected by buildings, objects and humans relating to absorption and reflection.

Kare Medikal ve Analitik Cihazlar Sanayi ve Ticaret LTD. ŞTİ.
Alıcı Osb Mh. 2017 cd. No:24
Sincan
Ankara, Turkey
Tel: +90 312 911 4 911
www.karemedical.eu

Imported and distributed by:
DeVilbiss Healthcare GmbH
Kamenzer Str. 3
68309 Mannheim
Germany

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42

Diagram of the Hikoneb 3000 ultrasonic nebulizer, showing the main unit and its components, numbered 1 through 38.

- 1: Base
- 3: Casters
- 6: Main column
- 8: Locking lever
- 9: Locking pin
- 10: Locking pin
- 11: Nebulizer chamber
- 12: Filter
- 13: Filter cap
- 14: Flexible delivery tube
- 15: Bypass valve
- 16: Bypass valve handle
- 17: Bypass valve handle
- 18: Bypass valve handle
- 19: Bypass valve handle
- 20: Bypass valve handle
- 21: Control unit
- 22: Filter cap
- 23: Filter cap
- 24: Filter cap
- 25: Filter cap
- 26: Filter cap
- 27: Filter cap
- 28: Filter cap
- 29: Filter cap
- 30: Filter cap
- 31: Filter cap
- 32: Filter cap
- 33: Filter cap
- 34: Filter cap
- 35: Filter cap
- 36: Filter cap
- 37: Filter cap
- 38: Filter cap

Inset detail of the control unit (21) showing the power button (4) and the display (5).

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