Manufacturer

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Monopolar electrodes

IMPORTANT!

Please read all information carefully.

These instructions for use do not replace the user manual of the electro-surgical unit used! Read the user manual of the electro-surgical unit and ask Erbe or your distributor in case of doubt!

1 Intended use

The monopolar electrodes are intended for monopolar cutting and/or coagulation of biological tissue.

2 Normal use

Number Group 21191-0xx: This monopolar electrode can be connected to electro-surgical units using electrode handles (shaft diameter: 4 mm).

Number Group 21191-1xx: This monopolar electrode can be connected to electro-surgical units using electrode handles (shaft diameter: 2.4 mm).

Erbe recommends the use of Erbe electrode handles. If a handle from another manufacture is used to operate the HF electrode, please observe the following: The joint from HF electrodes to handles of other manufacturers is not standardized. The handle of another manufacture, although it may appear to fit the Erbe HF electrode, may have different dimensions and tolerances, a different type of insulation, etc. Such differences hold risks for patient and staff. Erbe Elektromedizin is glad to provide more information as required.

Not intended for direct contact with the heart or the central circulatory system!

3 Maximum electrical capacity

The maximum electrical capacity of these products is:

- 4300 V_p

4 Safety instructions

WARNING!

This product may be used only by trained medical staff, who have been shown how to use it according to the instructions for use.

Check the product including all the insulation for mechanical damage and perform a function test every time before it is used.

If damaged, do not use this product!

Do not use in the presence of combustible or explosive materials!

Never lay this product on the patient or in his/her direct vicinity.

Protect the product from all physical damage! Do not throw! Do not apply force of any kind!

If components or parts of the product are bent, do not attempt to straighten them!

As an exception to this, the following products can be bent: item. nos. 21191-003, 21191-004, 21191-005, 21191-006, 21191-007, 21191-103, 21191-104, 21191-105, 21191-106, and 21191-107.

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

5 How to use

5.1 Only use cleaned, disinfected and sterilized products

Clean, disinfect and sterilize the product before using it for the first time and before every subsequent use.

5.2 Checking the product

Before each use, check the product including all the insulation for signs of physical damage!

If damaged, do not use this product!

5.3 Assembling the product

Insert the HF electrode into the electrode handle as it goes.

Before every use conduct a performance test.

6 Cleaning, disinfection, sterilization

6.1 Safety instructions

May be cleaned in an ultrasonic bath.

Never use sharp objects for cleaning.

Disinfectants must be rinsed off thoroughly after use.

Maximum air pressure during drying: 2 bar.

Maximum water pressure during rinsing: 2 bar.

During machine cleaning and disinfection, do not exceed 95 °C.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter "Overview of validated procedures".

During sterilization, do not exceed 138 °C.

Do not sterilize in hot air.

6.2 Processing limitation

Frequent processing can damage this product. Observe the safety instructions for checking the product. With visible damage or functional impairment, the product may no longer be used.

6.3 Dismantling

1. Remove the HF electrode from the handle.

6.4 Required aids

Processing step	Items Needed
Precleaning	Soft plastic brush/soft cloth
Manual cleaning/ disinfection	Soft plastic brush/soft disposable cloth (low particle) Compressed air gun (2 bar maximum)

6.5 Precleaning

Use water for precleaning, a non-fixing disinfectant if necessary.

Suitable for cleaning in a 3% hydrogen peroxide bath.

1. Remove surface contaminants with a soft brush/cloth while holding the product in a water bath and/or rinsing the product under running water.

6.6 Manual cleaning and disinfection

Use a liquid detergent suitable for the preparation of a dipping bath.

Use a disinfectant compatible with the detergent suitable for the preparation of a dipping bath.

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter "Overview of validated procedures".

1. Prepare a cleaning bath and a separate disinfectant bath.
2. Immerse the product completely in a cleaning bath, ensuring that the product does not come into contact with other components in the bath.
3. Thoroughly clean the surfaces with a soft brush/soft disposable cloth.
4. Rinse the surfaces with sterile, deionized water (at least 1 minute).
5. Repeat the preceding cleaning steps if the product still shows visible contamination.
6. Immerse the product completely in a disinfectant bath, ensuring that the product does not come into contact with other components in the bath.
7. Rinse the surfaces with sterile, deionized water (at least 1 minute).
8. Dry the product until there are no more fluid residues visible (e.g. with filtered compressed air).

6.7 Cleaning and disinfection by machine

The washer /disinfector must have fundamentally certified effectiveness (e.g. CE marking according to DIN EN ISO 15883).

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter "Overview of validated procedures".

1. Place the product in a suitable rinsing basket carefully. In doing so make sure the product is not touching other instruments or parts of instruments.

2. Start a tested program with the following properties:
 - Sufficient cleaning effect.
 - Thermal disinfection: 5 to 10 minutes at 90 to 95 °C, A₀ ≥ 3000.
 - Final rinse with distilled or deionized water.
 - Sufficient product drying (no visible fluid residues).
3. If contaminants are still visible on the product after the end of the routine, repeat precleaning and machine cleaning and disinfection.

6.8 Check

1. Check the product for visible damage and wear:
 - Damage to the product, e.g. cracks, rough surface, splintering.
 - Damage to the product insulation and/or cable/connector insulation, e.g. cracks and breaks.
 - Damage to the coating.
 - Bending of the product or parts of the product. Do not attempt to straighten bent instruments!

If damaged, do not use this product!

6.9 Packaging

1. Pack the product in disposable sterilization packaging (single or double packaging) made from paper/foil and/or in a sterilization container.

6.10 Sterilization

Only sterilize products that have been cleaned and disinfected.

Erbe Elektromedizin recommends steam sterilization with the method described below. If other sterilization methods are used, Erbe Elektromedizin shall bear no responsibility.

Steam sterilization

- Fractionated vacuum method with adequate product drying
- Holding time 3 to 18 minutes at 132 to 138 °C
- Sterilizers in accordance with applicable national standards and regulations (e.g. DIN EN 13060 or DIN EN 285)
- Sterilization process validated according to DIN EN ISO 17665

Please follow the recommendations issued by the sterilizer manufacturer concerning loading, handling, and drying times.

6.11 Overview of validated procedures

The following procedures have been validated as suitable for processing the product:

- Manual cleaning/disinfection with the detergent Cidezyme and disinfectant Cidex OPA (Johnson & Johnson Medical Limited, Gargrave, Skipton).
 - Cidezyme: Dosage 1.6%, 5 minutes exposure time at room temperature.
 - Cidex OPA: Dosage 0.3%, 12 minutes exposure time at 20 °C.
- Machine cleaning and disinfection in a G 7836 CD disinfector (Miele & Cie. KG, Gütersloh, Germany) with the cleaning agent needisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg, Germany).
 - needisher mediclean forte: Dosage 0.5%, 5 minutes exposure time at 55 °C
 - Disinfection: 5 minutes at 90 °C
- Sterilization in an autoclave Systec V-150 (Systec GmbH, Linden/Germany).
 - Steam sterilization with saturated steam, fractionated vacuum procedure, 3 minutes, 132 °C, 3 bar.

Erbe recommends the processing procedures described in these instructions for use. Equivalent different procedures are possible if not explicitly excluded. It is incumbent on the user to ensure the suitability of the actual procedures used by means of suitable measures (e.g. validation, routine monitoring, check of material compatibility).

7 Disposal

Dispose of the product, packaging material and accessories (if available) in accordance with the national guidelines and legislation applicable in each case.

8 Symbols

Symbol	Explanation
	Consult instructions for use
	Caution, consult accompanying documents

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© Erbe Elektromedizin GmbH