Manufacturer

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Bipolar instruments

IMPORTANT!

Please read all information carefully.

These instructions for use do not replace the user manual of the electrosurgical unit used! Read the user manual of the electrosurgical unit and ask Erbe or your distributor in case of doubt!

1 Intended use

The bipolar instruments are intended for bipolar coagulation.

2 Normal use

The bipolar instruments are connected to automatically controlled electrosurgical units.

3 Maximum electrical capacity

The maximum electrical capacity is given on the instrument in [V_p].

4 Safety instructions

WARNING!

This product may be used only by trained medical staff, who have been shown how to use it according to the instructions for use.

Before each use, check the product including all the insulation for signs of physical damage!

If damaged, do not use this product!

Do not use in the presence of combustible or explosive materials!

Always store the connected product isolated from the patient and user. Otherwise there is the risk of undesired coagulation, especially in the use of an AUTO START function.

Sharp instruments can cause injuries.

Protect this product from any form of mechanical damage! Do not throw! Do not use force!
The instrument tips heat up if subjected to heavy loads. In order to prevent unintentional coagulation, do not touch tissue with the instrument tips accidentally!

Erbe recommends using a smoke evacuator, protective goggles, and a respirator to protect against smoke.

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

5 How to use

Clean, disinfect and sterilize the product before using it for the first time and before every subsequent use.

5.1 Connecting the product

1. Attach the product to a matching Erbe connecting cable for bipolar instruments.
2. Connect the product to the bipolar socket of the electrosurgical unit.

5.2 Unit settings

Always work with the lowest power and voltage settings for the required surgical effect.

Only select unit settings that do not exceed the maximum electrical capacity of the product (V_p max or K_{Vp} max) as this could damage the product. Information on selecting appropriate unit settings is provided in the User Manual of the electrosurgical unit (for Erbe electrosurgical units, read the “Accessories” section).

Erbe electrosurgical unit: Erbe recommends operating the instrument in the soft coagulation mode. Sparks generated in the forced coagulation mode may cause increased wear.

5.3 During use

Regularly wipe blood and remnants of tissue from the instrument tips.

5.4 Precleaning the product in the procedure room

1. Use a soft cloth (such as a compress) to remove any obvious contamination immediately after use and/or place the product in deionized water.

6 Cleaning, disinfection, sterilization

6.1 Safety instructions

Suitable for cleaning in an ultrasonic bath (max. 10 minutes per processing cycle).

On no account use sharp or abrasive objects for cleaning purposes.

Disinfectants must be rinsed off thoroughly after use.

Maximum air pressure during drying: 2 bar.

Maximum water pressure during rinsing: 2 bar.

During machine cleaning and disinfection, do not exceed 95 °C.

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter “Overview of validated procedures”.

During sterilization, do not exceed 138 °C.

Do not sterilize in hot air.

Bipolar forceps: Do not sterilize with STERRAD (J&J)!

6.2 Reprocessing limitation

Frequent reprocessing and operative stress have impacts on this product. The service life of the product may be unexpectedly reduced depending on processing and type of application. Observe the safety instructions to inspect the product. No longer use a product that is obviously damaged or functionally impaired.

A minimum of 20 processing cycles should be assumed for normal use.

6.3 Dismantling

1. Disconnect the connecting cable from the product.

6.4 Required aids

Reprocessing step	Items Needed
Precleaning	Soft plastic brush/soft cloth
Manual cleaning/ disinfection	Soft plastic brush/soft disposable cloth (low particle) Compressed air gun (2 bar maximum)

6.5 Precleaning

Use water for precleaning, a non-fixing disinfectant if necessary.

Suitable for cleaning in a 3% hydrogen peroxide bath (max. 30 minutes per processing cycle).

1. Remove surface contaminants with a soft brush/cloth while holding the product in a water bath and/or rinsing the product under running water.

6.6 Manual cleaning and disinfection

Use a liquid detergent suitable for the preparation of a dipping bath.

Use a disinfectant compatible with the detergent suitable for the preparation of a dipping bath.

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3 (20195-003 and 20195-025; 5.5 to 11).

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter “Overview of validated procedures”.

1. Prepare a cleaning bath and a separate disinfectant bath.
2. Immerse the product completely in a cleaning bath, ensuring that the product does not come into contact with other components in the bath.
3. Thoroughly clean the surfaces with a soft brush/soft disposable cloth.
4. Rinse the surfaces with sterile, deionized water (at least 1 minute).
5. Repeat the preceding cleaning steps if the product still shows visible contamination.
6. Immerse the product completely in a disinfectant bath, ensuring that the product does not come into contact with other components in the bath.
7. Rinse the surfaces with sterile, deionized water (at least 1 minute).
8. Dry the product until there are no more fluid residues visible (e.g. with filtered compressed air).

6.7 Cleaning and disinfection by machine

The washer/disinfecter must have fundamentally certified effectiveness (e.g. CE marking according to DIN EN ISO 15883).

The detergent and disinfectant must be suitable for medical devices made of plastic and metal and have a pH value between 5.5 and 12.3 (20195-003 and 20195-025; 5.5 to 11).

Observe the specifications from the manufacturer of the detergent and disinfectant and in the chapter “Overview of validated procedures”.

1. Place the product in a suitable rinsing basket carefully. In doing so make sure the product is not touching other instruments or parts of instruments.
2. Start a tested program with the following properties:
 - Sufficient cleaning effect.
 - Thermal disinfection: 5 to 10 minutes at 90 to 95 °C, A₀ ≥ 3000.
 - Final rinse with distilled or deionized water.
 - Sufficient product drying (no visible fluid residues).
3. If contaminants are still visible on the product after the end of the routine, repeat precleaning and machine cleaning and disinfection.

6.8 Check

1. Check the product for visible contamination, visible damage and wear:
 - damage to the product, e.g. cracks, rough surface, splitting, corrosion.
 - Damage to the product insulation and/or cable/connector insulation, e.g. cracks and breaks.
 - Bending of the product or parts of the product. Do not attempt to straighten bent instruments!

If damaged, do not use this product!

Discolorations caused by corrosion can be removed by thoroughly cleaning with cleaning agents containing phosphoric acid. Only the tip of the instrument must come into contact with the cleaning solution.

6.9 Packaging

1. Pack the product in disposable sterilization packaging (single or double packaging) made from paper/foil and/or in a sterilization container.

6.10 Sterilization

Only sterilize products that have been cleaned and disinfected.

Erbe Elektromedizin recommends steam sterilization with the method described below. If other sterilization methods are used, Erbe Elektromedizin shall bear no responsibility.

Steam sterilization

- Fractionated vacuum method with adequate product drying
- Holding time 3 to 18 minutes at 132 to 138 °C
- Sterilizers in accordance with applicable national standards and regulations (e.g. DIN EN 13060 or DIN EN 285)
- Sterilization process validated according to DIN EN ISO 17665

Please follow the recommendations issued by the sterilizer manufacturer concerning loading, handling, and drying times.

6.11 Overview of validated procedures

The following procedures have been validated as suitable for processing the product:

- Manual cleaning/disinfection with the detergent Cidezyme and disinfectant Cidex OPA (Johnson & Johnson Medical Limited, Gargrave, Skipton).
- Cidezyme: Dosage 1.6%, 5 minutes exposure time at room temperature.
- Cidex OPA: Dosage 0.3%, 12 minutes exposure time at 20 °C.
- Machine cleaning and disinfection in a G 7836 CD disinfecter (Miele & Cie. KG, Gütersloh, Germany) with the cleaning agent neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg, Germany).
- neodisher mediclean forte: Dosage 0.5%, 5 minutes exposure time at 55 °C.
- Disinfection: 5 minutes at 90 °C.
- Sterilization in an autoclave Systec V-150 (Systec GmbH, Linden/Deutschland).
- Steam sterilization with saturated steam, fractionated vacuum procedure, 3 minutes, 132 °C, 3 bar.

Erbe recommends the processing procedures described in

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© Erbe Elektromedizin GmbH