MICRONAUT-S MRSA / GP

Nur zur *in vitro* Diagnostik

Deutsch

Revision A (Oktober 2018)

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TEST PRINCIPLE

The susceptibility testing with the MICRONAUT-S MRSA / GP test plate is based on the rehydration of antibiotics by adding a standardized bacterial suspension. After incubation for 18 - 24 hours at 35 - 37 °C the result is read with the plate reader and evaluated with the MICRONAUT Software or read visually and interpreted.

REAGENTS

Contents

- 1 test can be performed per each plate
- 40 test plates per packaging unit
- Sealing foils (not perforated)

Required additional reagents and materials:

- Mueller Hinton Broth, cation adjusted (CAMHB)
- MICRONAUT-H Medium
- Lysed horse blood for *Streptococcus pneumoniae* and other fastidious bacteria; Company Thermo Fisher Scientific (Oxoid) or Fiebig-Nährstofftechnik
- NaCl 0,9%
- 1-channel-reservoirs
- 8-channel-pipette (100-1200 µl) incl. pipette tips
- Plate reader (validated for MICRONAUT systems)
- MICRONAUT Software

Article numbers of reagents and materials are available at the local distributors.

Laboratory materials

- Incubator
- McFarland standard 0.5
- Blood agar plate (without additives)
- Inoculation loops
- Marking pen
- Pipette (50-1000 µl)

COMPOSITION OF THE MEDIA

Media	Components
NaCl 0.9%, 1l	Sodium chloride
Mueller Hinton Broth (cation-adjusted) 100 tubes of 11.5 ml (± 0.5 ml) 20 tubes of 11.5 ml (± 0.5 ml)	Beef Extract Acid Hydrolysate of Casein Starch
MICRONAUT-H Medium ¹ 100 tubes of 11.5 ml (± 0.5 ml)	Glucose Hematin Neopeptone Tween 80 Pyridoxal NAD Columbia Broth Base Yeast Extract Agarose Type II A

¹Compound confidential

Please note: The broth is chosen depending on the testing in accordance with international standards respective the bacteria group. Mueller Hinton Broth should be used as standard medium for testing non-fastidious Gram-positive bacteria. For particular bacteria groups it is advisable to add supplement to the broth.

Cation-adjusted Mueller-Hinton Broth has to be used for susceptibility testing of **daptomycin**.

STABILITY/ STORAGE

MICRONAUT-S MRSA / GP test plates have to be stored in the original packaging at 15 - 25 °C and can be used up to the indicated expiration date. Mueller Hinton Broth has to be stored at 15 - 25 °C, MICRONAUT-H Medium at 2 - 8 °C. For shelf life of the media please follow the expiration date printed on the product. For reagents follow the storage instructions indicated on the products.

PRECAUTIONARY MEASURES

- Only to be used for *in vitro* diagnostic.
- Do not pipette the reagents by mouth.
- Only for proper use.
- Samples, bacterial cultures and the inoculated test plates have to be considered as potentially infectious and must be treated properly and with respect to the corresponding precautionary measures by qualified specialist staff. It is important to work aseptically during the whole test procedure. For information please refer to "BioSafety in Microbiological and Biomedical Laboratories, HHS Publication No. (CDC) 99-8395, 4th Edition (April 1999)", or to the corresponding national legal requirements.
- A pouch filled with the indicator "silica gel" is supplied with the sealed test plates. The drying agent contains cobalt chloride. Please make sure not to damage the drying agent pouch. A colour change of the indicator pouch from blue to pink may indicate trapped moisture. Please do not use this plate.
- Upon reading and evaluation of the tests, all samples, inoculated and contaminated products (pipette tips, 1-channel-reservoir and test plates) must be autoclaved, burnt or disinfected in a bactericidal solution before disposal.
- It is important to follow the instructions carefully, each deviation may influence the quality of the results.
- The test results should be interpreted by qualified staff with experience in microbiology. The clinical background, origin of the samples, colony and microscopic morphology, serology and the identification result must be taken into consideration when interpreting the results.

TEST PROCEDURE

Preparation of the samples

- Prepare a tube with 1-5 ml NaCl 0.9%, pH 5.5 to 6.5.
- Prepare a tube of the corresponding Mueller Hinton Broth or MICRONAUT-H Medium.

Preparation of the inoculum

- Pick several single colonies of a maximum 24 hours old pure culture from the blood agar (without additives).
- Homogenize the colonies well in 1-5 ml NaCl 0.9 % until the turbidity matches a McFarland standard of 0.5.

Mueller Hinton Broth

- Gram-positive bacteria:
Pipette 100 µl of the bacterial suspension into 11.5 ml Mueller Hinton Broth and homogenize well.
- Fastidious bacteria (e. g. *Streptococcus pneumoniae* and other fastidious streptococci):
Pipette 200 µl of the bacterial suspension into 11.5 ml Mueller Hinton Broth. To ensure sufficient growth of fastidious bacteria lysed horse blood is added to the Mueller Hinton Broth. The final concentration of the lysed horse blood should be at least 2.5 % (275 µl to 11.5 ml Mueller Hinton Broth). Homogenize the supplemented medium well. The MICRONAUT-S MRSA / GP test plate must be read visually.
Remark: Photometric reading with MICRONAUT software **panel type „H“** is possible if exclusively lysed horse blood from Thermofisher Scientific (Oxoid) or Fiebig-Nährstofftechnik in the final concentration of 2.5% is used.

MICRONAUT-H Medium

- MICRONAUT-H Medium for fastidious bacteria (e. g. Coryneform bacteria, *Listeria*) and as an alternative test medium for streptococci:
Pipette 200 µl of the bacterial suspension into 11.5 ml MICRONAUT-H Medium and homogenize well.
Note: MICRONAUT-H Medium as test medium must not be used for MIC determination of daptomycin when testing fastidious bacteria (e. g. streptococci). Susceptibility testing of daptomycin requires the use of Mueller Hinton Broth.

Inoculation

- Remove the MICRONAUT-S MRSA / GP test plate from the packaging not more than 30 minutes before inoculation and dispose the drying agent.
- Label the test plate.
- Pour the prepared suspension into a 1-channel-reservoir.
- Inoculate the microtitration plate with 100 µl per well manually by using an 8-channel-pipette.

Sealing and incubation

Mueller Hinton Broth

- After the inoculation cover the test plates with the unperforated sealing foil or an unused microtitration plate (do not stack up more than 5 test plates).
- Place the test plate in an incubator at 35 - 37 °C for 18 - 24 hours.

MICRONAUT-H Medium

- After the inoculation stack up the test plates and cover the one on top with a microtitration plate or another lid. Do not cover with sealing foil (do not stack up more than 5 test plates).
- Place the test plate in an incubator at 35 - 37 °C for 22 - 24 hours in an enriched CO₂-atmosphere if necessary.

Reading

- Remove the sealing foils.
- Wipe off the bottom of the test plates.
- Read the test plates with the photometer or visually.

Evaluation

The growth control must be covered (turbid) otherwise the test must be repeated.

The values that are determined by photometric reading are evaluated, interpreted and the plausibility is checked with the MICRONAUT Software (see instructions for MICRONAUT Software). Test results are displayed on the screen or on the printout. Remarks may appear below the test results indicating an insufficient growth, MRSA phenotype or notifiable bacteria etc.

Results read visually should be documented on a printed layout or an evaluation protocol.

INTERPRETATION OF THE RESULTS

Breakpoint Interpretation: Antibiotics which are tested in two concentrations.

visual reading		automated reading		result	
1	2	1	2		
no growth	no growth	-	-	S	susceptible
growth	no growth	+	-	I	intermediate
growth	growth	+	+	R	resistant

1 = low breakpoint concentration

2 = high breakpoint concentration

Explanations

To judge the therapeutic susceptibility of bacteria, please refer to MIC interpretation standards of the corresponding guidelines (CLSI, EUCAST).

Enterococci strains

Check the MICRONAUT-S MRSA / GP test plates visually when testing enterococci.

MIC value

Minimum inhibitory concentration (MIC) is defined as the lowest concentration of an antibiotic that will inhibit visible growth of a microorganism. Bacterial growth can be observed as turbidity of the test medium within or as cell sediment (nodule formation) at the bottom of the respective wells of the microtitration plates. Prerequisite for the determination of MIC values is a positive growth control.

Skipped wells

Occasionally skipped wells could occur within the microtitration plates. Skipped wells appear due to lack of growth in one or more wells within an antibiotic concentration series while the wells of the next higher / lower concentration of the same antibiotic show bacterial growth.

Reasons for the occurrence of skipped wells can be complex: manifestation of hetero-resistance, inhomogeneous inoculum or inoculation, inadequate filling of the microtitration plates, cross contamination of cavities by carrying over of antibiotics during pipetting, working with mixed cultures, drops of bacterial suspension at the edge of the cavities and therefore no contact with the antibiotic, etc..

For MICRONAUT-S MRSA / GP test plates a single skipped well within a concentration series can be ignored. The concentration of the well from which on no further growth occurs, can be read as the MIC value. If several skipped wells occur in a concentration series the test must be repeated. If skipped wells occur in breakpoint concentrations (maximum of two concentrations per antibiotic agent) the test must be repeated.

EXPLANATION OF THE ANTIBIOTIC CONFIGURATION OF THE MICRONAUT-S MRSA / GP TEST PLATE

The antibiotic configuration of the MICRONAUT-S MRSA / GP test plate allows the specific detection of clinically relevant single or multi-resistances of Gram-positive bacteria causing nosocomial infections. The susceptibility testing with highly effective reserve antibiotics offers alternatives in case of extreme multiresistance.

Susceptibility testing of staphylococci

- Detection of staphylococcal penicillinases by MIC determination of penicillin G
- Detection of oxacillin resistance by MIC determination of oxacillin and ceftiofur
- Detection of inducible MLS_B resistance by testing the antibiotic combination erythromycin/clindamycin
- Detection of phenotypic resistance patterns of epidemic MRSA
- Testing of novel antibiotics

Susceptibility testing of enterococci

- Detection of ampicillin resistance by MIC determination of ampicillin
- Detection of phenotypic glycopeptide resistance patterns of vancomycin-resistant enterococci by determining the MIC of teicoplanin and vancomycin
- Detection of HLAR-strains through high-level-resistance testing of gentamicin
- Testing of novel antibiotics

Susceptibility testing of pneumococci

- Detection of PBP-changes by MIC determination of penicillin G
- Detection of group IV quinolones resistance by MIC determination of moxifloxacin
- Detection of erythromycin resistance
- Detection of vancomycin resistance
- Testing of novel antibiotics

TECHNICAL REMARKS

In order to obtain best results please follow the below listed points of the instructions carefully:

- MICRONAUT-S MRSA / GP test plates are intended for single use only. Test plates cannot be reused.
- Work with pure culture from blood agar (without additives) not older than 24 hours.
Exception: Fastidious bacteria can be taken from a blood agar plate (without additives) after 48 hours of incubation at 35 - 37 °C.
- If you do not use culture media ready for use follow the instruction carefully when preparing the broth.
- When using Mueller Hinton Broth, check the correct concentration of the bivalent cations Ca²⁺ and Mg²⁺ (20 to 25 mg Ca²⁺/l and 10 to 12.5 mg Mg²⁺/l). Verify the pH-value of the Mueller Hinton Broth. It should be between 7.2 and 7.4 at room temperature (25 °C).
- Use NaCl 0.9 %, pH 5.5 - 6.5.
- Store the MICRONAUT-H Medium at 2 – 8 °C.
- Preheat the medium before using it (1 hour in the incubator).
- Follow the correct McFarland standard 0.5 adjustment of the NaCl suspension. Homogenize the suspension sufficiently.
- Transfer the bacteria suspension on a blood agar plate (without additives) for purity control.
- Follow the incubation times carefully. Do not incubate non-fastidious bacteria less than 18 hours and fastidious bacteria less than 22 hours.

QUALITY CONTROL

The test plates and reagents are subject to quality controls which are carried out systematically at different stages of the production. The bacteriological quality control can be carried out with the following strains.

Strains	ATCC No. ¹	DSMZ No. ²
<i>Staphylococcus aureus</i>	ATCC 29213	DSM 2569
<i>Enterococcus faecalis</i>	ATCC 29212	DSM 2570
<i>Streptococcus pneumoniae</i>	ATCC 49619	DSM 11967

¹ATCC = American Type Culture Collection

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QUALITY AND PERFORMANCE DATA

The general requirements apply to DIN EN ISO 20776-1. For monitoring the accuracy for minimal inhibitory concentrations (MICs) please refer to the acceptable limits for quality control strains according to the relevant version of CLSI document M100.

Attention: Deviations from the acceptable QC MIC ranges can occur with **fosfomycin** antimicrobial susceptibility testing (AST) due to expression of hetero-resistance based of resistant subpopulations. These deviations are linked to the high inoculum size of the broth microdilution procedure (refer to Ballesterro-Tellez et al., 2017).

GUARANTEE

The quality data of the MICRONAUT-S MRSA / GP test plate has been determined by strictly following the present instruction. Divergences or alterations of the test procedure may reduce the quality of the results. Any claims for damages are excluded in this case.


EXPLANATION OF THE SYMBOLS ON THE LABELS

	For single use only
	Number of possible tests
	Storage conditions
	Instructions for use
	Security advice on the safety data sheet
	Expiry date
	CE marking in accordance with 98/79/EC (IVDD)
LOT	Indication of the lot number
IVD	In vitro diagnostics
REF	Article number

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SHORT INSTRUCTION MICRONAUT-S MRSA / GP

* MICRONAUT-H Medium for fastidious bacteria must not be used for MIC determination of daptomycin. Instead, cation-adjusted Mueller Hinton Broth (CAMHB) supplemented with 2.5% lysed horse blood must be used.

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