

Oxi/FermPluri - Test

REF 78620 / 78621

DATE		SAMPLE														
	GROUP 1			GROUP 2			GROUP 3			GROUP 4			GROUP 5			
Test	Anaerobic Glucose	Arginine	Lysine	Lactose	N ₂	Sucrose	Indole	Xylose	Aerobic Glucose	Maltose	Mannitol	Phenyl - alanine	Urea	Citrate	Oxidase	
Positivity code	4	2	1	4	2	1	4	2	1	4	2	1	4	2	1	
Results																
Code sum																
NUMERICAL CODE		MICROORGANISMS														

REACTIONS TABLE

Negative reactions ⊖

Anaerobic Glucose Arginine Lysine Lactose - N₂ Sucrose - Indole Xylose Aerobic Glucose Maltose Mannitol Phenyl - alanine Urea Citrate

Positive reactions ⊕

Anaerobic Glucose Arginine Lysine Lactose - N₂ Sucrose - Indole Xylose Aerobic Glucose Maltose Mannitol Phenyl - alanine Urea Citrate

N₂ Indole

Oxi/FermPluri - Test

REF 78620 / 78621

DATE		SAMPLE														
	GROUP 1			GROUP 2			GROUP 3			GROUP 4			GROUP 5			
Test	Anaerobic Glucose	Arginine	Lysine	Lactose	N ₂	Sucrose	Indole	Xylose	Aerobic Glucose	Maltose	Mannitol	Phenyl - alanine	Urea	Citrate	Oxidase	
Positivity code	4	2	1	4	2	1	4	2	1	4	2	1	4	2	1	
Results																
Code sum																
NUMERICAL CODE		MICROORGANISMS														