

Merlin Diagnostika GmbH
Kleinstraße 14
D – 53332 Bornheim

Certificate of Analysis

Date of issue: 2019-04-15
Article number: E2-338-100
Material description: MICRONAUT – ID Medium
Lot: 3280319

Date of manufacture: 2019-04-11
Expiry date: 2020-04-11
Batch status: Accepted
Temp. Cond.: + 15 °C to + 25 °C

Composition:

Sodium chloride

This medium can be adjusted and/or supplemented according to the performance criteria required.

Test	Specification	Result
Appearance	clear, colourless	Conforms
pH-value	4.5 – 8.0	6.98
Filling Volume	5.5 ml ± 0.5 ml	Conforms
Sterility control min. 72 h; 30-35 °C	no growth	Conforms

This document has been established electronically and is valid without signature.