




































	Enterics	Pseudo- monaden	Staphylo- kokken	Entero- kokken	Strepto- kokken ^z	Acineto- bacter	Stenotropho- monas
CAZ 10µg (Ceftazidim)							
Ring- größen	19 / 22	17 / 17 ^o	X	X	X	X	X
CIP 5µg (Ciprofloxacin)							
Ring- größen	24 / 26 ^D	26 / 26	24 / 24 ^{L, MM}	15 / 15	X	21 / 21	X
NOR 10µg (Norfloxacin)							
Ring- größen	19 / 22 ^C	X	17 / – ^L	12 / 12	12 / –	X	X
GM 10 µg (Gentamicin)							
Ring- größen	14 / 17	15 / 15	22 / 22 ^M	X	X	17 / 17	X
FOT 50 µg (Fosfomycin)							
Ring- größen	E	X	X	X	X	X	X
SXT 23,75/1,25µg (Cotrimoxazol)							
Ring- größen	13 / 16	X	14 / 17	21 / 50	15 / 18	13 / 16	16 / 16 ^A

Stand: 25.01.2017

Ref: EUCAST v 7.0 (01.01.17) / NAK



	Enterics	Pseudo- monaden	Staphylo- kokken	Entero- kokken	Strepto- kokken ^z	Acineto- bacter
SAM 10/10µg (Ampicillin / Sulbactam)						
Ring- größen	14 / 50 ^A	X	H, I	B	U	X
OX 1µg (Oxacillin)						
Ring- größen	X	X	H	X	X	X
CXM 30µg (Cefuroxim)						
Ring- größen	19 / 19 ^c	X	J	X	U	X
MEM 10µg (Meropenem)						
Ring- größen	16 / 22	18 / 24	J	X	U	15 / 21
LEV 5µg (Levofloxacin)						
Ring- größen	19 / 22	22 / 22	24 / 24 ^{L, MM}	15 / 15 ^L	17 / 17	20 / 23
NI 100µg (Nitrofurantoin)						
Ring- größen	11 / 11 ^{C, F}	X	13 / 13 ^{C, P}	15 / 15 ^{C, S}	15 / 15 ^{C, V}	X

















Stand: 25.01.2017

Ref: EUCAST v 7.0 (01.01.17) / NAK



	Enterics	Pseudo- monaden	Staphylo- kokken	Entero- kokken	Strepto- kokken ^z	Acineto- bacter	Stenotropho- monas
PG 1IU (Penicillin)							
Ring- größen	X	X	26 / 26 ^{6,H,T}	X	18 / 18	X	X
AP 2µg (Ampicillin)							
Ring- größen	14 / 50 ^{A, AA} Für Enterics ist 10µg vorgeschrieben.	X	18 / 18 ^I	8 / 10	U	X	X
CIP 5µg (Ciprofloxacin)							
Ring- größen	19 / 22 ^D	22 / 25	24 / 24 ^{L, MM}	15 / 15 ^Q	X	21 / 21	X
NOR 10µg (Norfloxacin)							
Ring- größen	19 / 22	X	17 / - ^L	12 / 12	12 / -	X	X
T 30µg (Tetracyclin)							
Ring- größen	X	X	19 / 22 ^O	X	20 / 23	X	X
SXT23,75/1,25µg (Cotrimoxazol)							
Ring- größen	13 / 16	X	14 / 17	21 / 50	15 / 18	13 / 16	16 / 16 ^A



	Enterics	Pseudo- monaden	Staphylo- kokken	Entero- kokken	Strepto- kokken ^z	Acineto- bacter
CPD 10µg (Cefpodoxim)						
Ring- größen	21 / 21 ^c	X	J	X	U	X
LEV 5µg (Levofloxacin)						
Ring- größen	19 / 23	17 / 20	24 / 24 ^{L, MM}	15 / 15	17 / 17	20 / 23
OFX 5µg (Ofloxacin)						
Ring- größen	22 / 24	X	24 / 24 ^{L, MM}	X	X	X
E 15µg (Erythromycin)						
Ring- größen	X	X	18 / 21	X	18 / 21	X
NI 100µg (Nitrofurantoin)						
Ring- größen	11 / 11 ^{c, F}	X	13 / 13 ^{c, P}	15 / 15 ^{c, S}	15 / 15 ^{c, V}	X
FOX 30µg (Cefoxitin)						
Ring- größen	19 / 19	X	22 / 22 ^{K, KK}	X	X	X

Stand: 25.01.2017

Ref: EUCAST v 7.0 (01.01.17) / NAK